


## **Surgical Instrument Trays 100% Stainless Steel Surgical Instrument Trays prevent rusting and corrosion.**

Handmade German craftsmanship renders high quality Surgical Instrument Trays. ergo-Logistics Surgical Instrument Trays are 100% high-tensile-strength electro-polished Stainless Steel. Each Surgical Instrument Tray is inspected by hand for smooth welding. ergo-Logistics Instrument Trays are designed to deliver optimal flow-thru for effective presoaking and cleaning of surgical instruments, Optimal flow through and venting of sides and bottom deliver the maximum penetration and venting for use with surgical instrument washers and steam sterilizers.

### **Surgical Instrument Trays – Wire Mesh rounded sides & ‘drop’ handles, all stainless steel construction**

<b>Product CODE</b>	<b>Description</b>	<b>L x W x H inches</b>
7500053	Surgical Instrument Tray Smooth wire mesh rounded corners	9.5 x 10 x 2
7500001	Surgical Instrument Tray Smooth wire mesh rounded corners	9.5 x 10 x 2.5
7500004	Surgical Instrument Tray Smooth wire mesh rounded corners	19 x 10 x 2
7500059	<u>Surgical Instrument Tray</u> Smooth wire mesh rounded corners	19 x 10 x 2.5
7500049	Surgical Instrument Tray Smooth wire mesh rounded corners	19 x 10 x 4
7500065	Surgical Instrument Tray Smooth wire mesh rounded corners	20 x 10 x 2
7500070	<u>Surgical Instrument Tray</u> Smooth wire mesh rounded corners	21 x 10 x 2.5
7501158	Surgical Instrument Tray Smooth wire mesh rounded corners	10 x 6.5 x 2.7
7501159	Surgical Instrument Tray Smooth wire mesh rounded corners	13.4 x 10 x 2.7
7501160	<u>Surgical Instrument Tray</u> Smooth wire mesh rounded corners	17.7 x 13.4 x 2.7




### **Surgical Instrument Trays – no handles, STACKABLE**

7500046	Instrument Tray Smooth wire mesh rounded corners	9.5 x 10 x 3
7500048	Instrument Tray Smooth wire mesh rounded corners	19 x 10 x 2

### **Surgical Instrument Trays – Perforated Steel Sides flat wire base for fast drying - with ‘drop handles’ sides are constructed from Stainless Steel and have .05” perforations every .23”**

<b>Product CODE</b>	<b>Description</b>	<b>L x W x H inches</b>
7500599	Surgical Instrument Tray Perforated sides flat wire bottom	19 x 10 x 2.5
7500729	Surgical Instrument Tray Perforated sides flat wire bottom	19 x 12 x 2.5
7515006	Surgical Instrument Tray Perforated sides flat wire bottom	19 x 10 x 4
7515007	<u>Surgical Instrument Tray</u> Perforated sides flat wire bottom	14 x 10 x 2.5
7500645	Surgical Instrument Tray Perforated sides flat wire bottom	9.5 x 10 x 2.5
7500809	Surgical Instrument Tray Perforated sides flat wire bottom	4.7 x 10 x 2.5
7500684	<u>Surgical Instrument Tray</u> Perforated sides flat wire bottom	18 x 12.5 x 2.5
7500730	Surgical Instrument Tray Perforated sides flat wire bottom	9 x 12.6 x 2.6
7500683	<u>Surgical Instrument Tray</u> Perforated sides flat wire bottom	4.5 x 12.6 x 2.5




### **Surgical Instrument Trays – Perforated Sides flat wire base for fast drying - with ‘combo handles’**

<b>Product CODE</b>	<b>Description</b>	<b>L x W x H inches</b>
7501136	Surgical Instrument Tray Perforated Sides	9.5 x 10 x 2.5
7501137	Surgical Instrument Tray Perforated Sides	19 x 10 x 2.5



### **Endoscope Baskets - twist locking lid and silicone "hold-downs", all stainless steel construction**

<b>Product CODE</b>	<b>Description</b>	<b>L x W x H inches</b>
7500849	Endoscope [1 scope] Basket with top locking latch	25 x 6 x 3
7500479	Endoscope [1 scope] Basket with top locking latch	11.5 x 3 x 2
7500475	Endoscope [1 scope] Basket with top locking latch	18 x 3 x 2
7500508	Endoscope [1 scope] Basket with top locking latch	26.5 x 3 x 2
7500831	Endoscope [1 scope] Basket with top locking latch	11.5 x 3 x 2
7500830	Endoscope [2 scopes] Basket with top locking latch	18 x 3 x 2



**Customer Service: Joan Hinchman, 1-800-341-7048, Fax: 1- 800 - 704 - 2815 email: [CustomerService@ergo-Logistics.com](mailto:CustomerService@ergo-Logistics.com),  
Technical Support: John Temple, 1- 509-747- 5027 [Support@ergo-Logistics.com](mailto:Support@ergo-Logistics.com)**